

User Guide

Medit Link



Medit Link **User Guide**

V2.3.0

Updated: 2020.09

Contents

1.	Introduction and Overview	7
1.1	Cloud Service	9
1.2	Workflow management using ordering system.....	9
1.3	Utilizing the Dashboard	9
1.4	Installation.....	10
1.4.1	System Requirements.....	10
1.4.2	Installation Guide.....	10
2	Registration	13
2.1	Landing Page	13
2.2	Registration.....	13
2.3	Login and Welcome Page	14
3	Medit Link Web	15
3.1	Dashboard	15
3.2	Communicating via Case Talk	16
3.2.1	Case Talk.....	16
3.2.2	Creating a Case Talk.....	17
3.2.3	Web Viewer.....	19
3.2.4	Model Tree.....	22
3.3	Partners.....	24
3.3.1	Partner Search	24
3.3.2	Recommended Partners.....	24
3.3.3	My Partner List	25
3.3.4	Partner Invitation	25
3.4	Membership	26
3.5	Redeem Code	26
3.6	Payment.....	26
3.6.1	Payment Details (Clinics, Labs accounts)	26

3.6.2	Payment account information (Lab account only)	27
3.7	My Information	27
3.7.1	My Account	27
3.7.2	My Organization (Admin account only).....	28
3.7.3	My Price List (Lab: Admin and Member accounts)	28
3.8	Help Centre	29
3.9	Announcements	29
3.10	Notices and Updates	29
3.11	Download the App	29
4	Medit Link App Overview	30
4.1	Title Bar	30
4.1.1	Network Status	30
4.1.2	Cloud Synchronization.....	31
4.1.3	User Information.....	32
4.1.4	Notifications.....	33
4.2	Boxes Overview	34
4.3	Settings	35
4.3.1	Patient Management.....	40
4.3.2	Notices and Updates.....	40
4.3.3	My Account	40
4.3.4	Help Center	40
4.4	Detailed Functions of Case Box/Work Box	41
4.4.1	Creating cases	41
4.4.2	Change viewing options (Clinic Account only)	41
5.	Medit Link App – Box Detail	43
5.1	Case Box / Work Box	43
5.2	Order Box / In Box	45
5.3	App Box	46
5.4	Trash Box.....	47

6. Detailed Explanation of the Process.....	48
6.1 Managing cases: creating, editing, cloning, deleting	48
6.1.1 Creating cases	48
6.1.2 Editing cases	50
6.1.3 Cloning cases.....	51
6.1.4 Deleting cases (Lab and Clinic).....	52
6.2 Inputting form information	53
6.2.1 Inputting form information in Case Box	53
6.2.2 Inputting form information in Work Box	53
6.2.3 Inputting bridge information	55
6.2.4 Setting scanning options.....	55
6.3 Copying form information.....	56
6.4 Deleting treatment information.....	57
6.5 Adding product information to 'Favorite'	58
6.5.1 Adding to favorites	58
6.5.2 Managing favorite templates.....	58
6.6 Copying form information between Scan and CAD tabs	59
6.6.1 From Scan tab to CAD tab.....	59
6.6.2 From CAD tab to Scan tab.....	59
6.7 Scanning process.....	60
6.8 CAD process.....	61
6.9 Exporting result files / Case Talk.....	62
6.9.1 Export.....	62
6.9.2 Case Talk.....	63
6.10 File Viewer.....	64
6.11 Model Tree	64
6.11.1 Covering empty spaces	66
6.11.2 Capture and Memo	66
6.11.3 Orthodontic Capture	67

6.12	Order process and order management	67
6.12.1	Ordering.....	67
6.12.2	Job Completed – Work Box.....	70
6.12.3	Checking order information	72
6.12.4	Managing expected/requested delivery dates	72
6.13	Patient Management	75
6.13.1	Adding patients.....	75
6.13.2	Deleting patients.....	76
6.13.3	Merging patients	76
6.14	Case filtering.....	77
6.14.1	Filtering by status	77
6.14.2	Filtering using period settings	78
6.14.3	Groping cases.....	79
6.14.4	Filtering with search capabilities.....	80

1. Introduction and Overview

Medit Link is a collaborative tool designed to facilitate communication between clinics and labs to support patient-centered workflows. It includes tools according to your business type (Clinic or Lab) and specific features can be used according to the account type. Medit Link helps users manage a series of digital dentistry workflows including scanning with Medit's intraoral and table top scanners, prosthesis order and manufacturing process for prosthesis through a user-friendly system.

Medit Link consists of a web-service and a computer application.

Medit Link Overview

		Web	App
		Register, search and manage partners, share cases via Web Viewer with a sharable link; access cases made through the App.	Scan, send/receive orders to/from the partners found through Web, enter form information; do CAD and manage patients.
Clinic	Admin	Search and manage partners, as well as your organization. <ul style="list-style-type: none"> • Dashboard • Case Talk • Partners - Search for Partners, Featured Partners, My Partner List • Membership • Redeem Code • Payment • My Info - My Account Info • My Info - My Organization Info. (Address, Organization, Additional, Scanner) • Help Center 	Scan, do CAD, send orders. <ul style="list-style-type: none"> • Dashboard • Case Box • Order Box • Trash Box - Permanently Deleting, Restore • App Box • Patient Management

	Member	<ul style="list-style-type: none"> • Case Talk • My Info - My Account • Help Center 	<ul style="list-style-type: none"> • Case Box • Trash Box - Restore App Box
Lab	Admin	<p>Search and manage partners, as well as your organization.</p> <ul style="list-style-type: none"> • Dashboard • Case Talk • Partners - Search For Partners, Featured Partners, My Partner List • Membership • Redeem Code • Payment • My Info - My Account <ul style="list-style-type: none"> ◦ My Info - My Organization ◦ My Info - My Price List • Help Center 	<p>Scan, do CAD, receive orders.</p> <ul style="list-style-type: none"> • Dashboard • In Box • Work Box • Trash Box - Permanently Deleting, Restore • App Box
	Member	<p>Scan and share cases.</p> <ul style="list-style-type: none"> • Case Talk • My Info - My Account • Help Center 	<p>Scan and do CAD.</p> <ul style="list-style-type: none"> • Work Box • Trash Box - Restore • App Box

1.1 Cloud Service

Medit Link supports automatic upload/download of database and files through AWS (Amazon Web Services). The data is automatically synchronized with computers of multiple users in the organization so that they can access this data from any computer at any time. Being connected to the Cloud helps with data backup.

The Cloud is also useful when it comes to quickly and easily sharing data between the dentist and the labs, creating a working environment that minimizes the hassle of sharing data using email. All data patient-related information is securely encrypted, stored and transmitted, and complies with medical law, HIPAA regulations, and GDPR regulations.

1.2 Workflow management using ordering system

Medit Link provides order and payment system.


You can quickly and easily share your data by signing a contract between the clinic and the lab.

You can view/manage the status of your orders.

Credit card-backed system works in conjunction with online payment platform companies Stripe (www.stripe.com) and I'mPort (<http://www.iampor.kr>) payment modules, making payment services easy and convenient.

1.3 Utilizing the Dashboard

You can view the statistics for the entire work done by Medit Link in a visualized graphs and chart – both in App and Web versions of Medit Link.



1.4 Installation


1.4.1 System Requirements

CPU	Intel Core i5 2.6GHz or higher
RAM	16 GB or higher
Graphic	NVIDIA GeForce GTX 760 (2GB) or higher / or equivalent AMD video card
OS	Window 8 64 Bit (unavailable in 32 Bit) or higher

1.4.2 Installation Guide

Medit Link is installed as a package with Medit Scan for Clinics and Medit Scan for Labs.

- Run Medit_Link_X.X.X.X.exe.



Medit Link User Guide

- Select the language for installation.




- Select the folder for installation of Medit Link. Read "License terms & conditions". carefully and Check the box "I agree to the License terms & Conditions." Click "Install" to proceed to the installation.



- As Medit Link is installed with Medit Scan for Clinics, the installation will not proceed if i500 device is connected to the computer.

- Make sure that the scanner cable is disconnected from the computer before proceeding to the next step.



- From this step, the installation will continue automatically and may take a few minutes. Do not turn off the computer until the installation is complete.



- After the installation is complete, restart the computer to use Medit Scan for Clinics and Medit Scan for Labs with Medit Link. You can skip this step If you do not have a scanner.




Note: the scanners may not work without restarting the computer.

2 Registration

2.1 Landing Page

Register at www.meditlink.com homepage.



2.2 Registration

Clinic	Lab
Administrator Dentist (Manager) All Medit Link functions are available. <ul style="list-style-type: none">✓ Order / Payment✓ Member Management Register	Administrator Representative / Leader (Manager) All Medit Link functions are available. <ul style="list-style-type: none">✓ Price Management for Dental Prostheses✓ Order Management✓ Member Management Register
Member Dental Hygienist and Dental Staff You can join your clinic's group as a member to use Medit Link. Case management and basic functions are available. <ul style="list-style-type: none">✓ Oral Scanning✓ Ordering of Dental Prostheses Register	Member Laboratory Employees You can join your lab's group as a member to use Medit Link. Case management and basic functions are available. <ul style="list-style-type: none">✓ Receiving Data for Making Dental Prostheses✓ Sharing Work Files and Production Processes Register

Medit Link requires an initial registration as Clinic (administrator & member) or Lab (administrator & member).

Administrator is usually a representative of your organization and has unlimited access to all the services, whereas, members have a limited access to Medit Link services.


- After specifying the name of your clinic or lab, click **"Check Availability"**.
 - If you are registering as a member, your organization should be activated prior to your

registration.

- Once you complete all the fields and click 'Register', a verification email will be sent to you.
- Verify your email address through the link in the email to complete the registration process.

2.3 Login and Welcome Page

- On the login page, enter your email address and password and click login.
 - To use Medit Link services, information about your organization should be entered using the administrator account.
 - The information entered on the **Welcome Page** (My Organization etc.) can be edited anytime on 'My Organization' page.
 - You can also register and manage intraoral and table top scanners and milling machines.
 - If you agreed to receive marketing information during registration, you will receive information about the registered products.



3 Medit Link Web

Case Talk
You can share the cases created in Medit Link and communicate with the patient through Web Viewer.

2020-07-21 ~ 2020-08-20 Case Name Search Case Name

Status	Case Name	Patient Name	Date of Creation	Scanned Date	Date of Sharing	Form Information
	test's Case - clone - clone	test	2020-08-20 18:09:31			13,14,15,16 Case Talk
	test's Case - clone	test	2020-08-20 17:15:38			13,14,15,16 Case Talk
	te's Case	te	2020-08-18 10:06:25			Case Talk
	test11231	test	2020-08-18 10:02:39			11,12,13,14,15,16 Case Talk
	test's Case - clone	test	2020-08-17 16:08:31			Case Talk
	test's Case	test	2020-08-17 16:04:47			13,14,15,16 Case Talk
	te's Case	te	2020-08-17 12:54:34			Case Talk
	tesw's Case	tesw	2020-08-14 18:41:18			Case Talk
	test's Case	test	2020-08-14 18:40:53			Case Talk
	test's Case	test	2020-08-14 18:34:39			25 Case Talk

Download Medit Link

3.1 Dashboard

Scanner
I500 (AA1803100019)
Usage Time: 11 hours 38 minutes
Calibration Period: 14 days
Elapsed Days: 371 days

Cloud Storage
TOTAL 5TB
Space Used: 0.86 GB
Available Space: 5,119.14 GB

Statistics
Total Cases: 11
Scanned Cases: 1
Scanning Time: 00:00:00
Total Orders: 2
Completed Orders: 0
Canceled Orders: 1
Rejected Orders: 0
Purchase Amount (USD): 0.00

Order Status
Pending: 1
Accepted: 0
Shipped: 1
Completed: 0

Graph
2020-08-10 ~ 2020-09-08
Total Cases: 19
Scanned Cases: 0
Scanning Time: 00:00:00
Total Orders: 0
Purchase Amount: 0

Tooth Type
2020-08-10 ~ 2020-09-08
Inlay: 1

Graphs and charts show statistics for the entire work done in Medit Link.


Scanner	Provides information about the active scanners, including scanner usage time, calibration interval, and calibration elapsed time.
Cloud Storage	Displays the total capacity of the cloud and the amount in use.
Statistics	Provides statistics on cases based on job status.
Order Status	Provides statistics on cases based on order status.
Graph	Graph job information by period.
Tooth Type	Provides statistics on the number of scans and orders by treatment information.


3.2 Communicating via Case Talk

3.2.1 Case Talk

Cases created in the Medit Link application can be shared with without any restrictions, or with specific people, and registered members of the Medit Link.

► How to



 Shared cases can be viewed via a browser without installing a separate program, enabling easy communication over the 3D data on various devices.

Medit Link User Guide

The screenshot shows the Medit Link Case Talk interface. On the left is a sidebar with navigation links: DASHBOARD, CASE TALK (selected), PARTNERS, MEMBERSHIP, REDEEM CODE, PAYMENT, MY INFO, and HELP CENTER. Below these are links for '새소식' (New News), 'Notice & Update', and a 'Download Medit Link' button. The main area is titled 'Case Talk' with a subtitle 'You can share the cases created in Medit Link and communicate with the patient through Web Viewer.' It includes a date range filter (2020-07-21 ~ 2020-08-20), a 'Case Name' dropdown, and a search bar. A table lists cases with columns: Status, Case Name, Patient Name, Date of Creation, Scanned Date, Date of Sharing, and Form Information. Each row has a 'Case Talk' button. The table contains 11 rows of test cases.

Status	Case Name	Patient Name	Date of Creation	Scanned Date	Date of Sharing	Form Information	
	test's Case - clone - clone	test	2020-08-20 18:09:31			13,14,15,16	Case Talk
	test's Case - clone	test	2020-08-20 17:15:38			13,14,15,16	Case Talk
	te's Case	te	2020-08-18 10:06:25				Case Talk
	test11231	test	2020-08-18 10:02:39			11,12,13,14,15,16	Case Talk
	test's Case - clone	test	2020-08-17 16:08:31				Case Talk
	test's Case	test	2020-08-17 16:04:47			13,14,15,16	Case Talk
	te's Case	te	2020-08-17 12:54:34				Case Talk
	tesw's Case	tesw	2020-08-14 18:41:18				Case Talk
	test's Case	test	2020-08-14 18:40:53				Case Talk
	test's Case	test	2020-08-14 18:34:39			25	Case Talk

3.2.2 Creating a Case Talk

Create a shared link from a case created by the Medit Link application.

- Select the case you want to share.

The screenshot shows the Medit Link Case Talk interface with a different date range (2020-08-11 ~ 2020-09-09). The table lists cases with columns: Status, Case Name, Patient Name, Form Information, Date of Creation, Scanned Date, Date of Sharing, and a 'Share' button. The first row, 'test's Case - clone' with patient 'testq1231', is highlighted with a red border. The table contains 11 rows of test cases.

Status	Case Name	Patient Name	Form Information	Date of Creation	Scanned Date	Date of Sharing	
	test's Case - clone	testq1231	22,23,24,25,26,27	2020-09-09 18:15:00			Share
	test's Case	testq1231	12,13,14,15,16, 23,24,25,26	2020-09-09 17:56:34			Share
	test's Case	testq1231	24,25,26	2020-09-09 17:53:24			Share
	test's Case	testq1231	25	2020-09-09 14:46:13			Share
	test's Case	testq1231	11,12,13,14,15,16,17,18,21,22,23,24	2020-09-08 18:44:48			Share
	test's Case	testq1231		2020-09-08 18:42:04			Share
	test's Case	testq1231	22,23,24,25,26,27	2020-09-08 17:42:05			Share
	test's Case	testq1231	24,25,26,27	2020-09-08 17:41:38			Share
	test's Case	testq1231	16	2020-09-08 17:38:51			Share
	test's Case	testq1231	11,12,13,14,15,16	2020-09-08 17:24:12			Share

- On the next screen, select the file you want to share.

Share to Case Talk

Case Talk

Case Talk List

Case Talk Count : 0

File List

Get Link

Name	Type	Size
<input checked="" type="checkbox"/> Maxilla	meditMesh	1,133 KB
<input checked="" type="checkbox"/> Default_1.stl	stl	2,694 KB
<input checked="" type="checkbox"/> Default_2.stl	stl	2,596 KB
<input checked="" type="checkbox"/> Default_3.stl	stl	2,593 KB
<input type="checkbox"/> Default_4.stl	stl	2,296 KB
<input type="checkbox"/> Default_5.stl	stl	2,430 KB
<input type="checkbox"/> Default_6.stl	stl	2,609 KB

4 items selected

8.81MB / 200MB

ⓘ If you share files over 200MB, Web Viewer may not work smoothly on some mobile devices.
 Please check the recommend requirement.
[Check requirement](#)

ⓘ You can create Case Talk and share the link with anyone, specific people, organization in Medit Link.
 You can set the password through website to limit the access to the shared link.
[Go to Case Talk](#)

⚠ Please note that Medit Link is not responsible for any problems caused by sharing the public link you created.

Close

- Click "Get Link" to create a shared link.

Share to Case Talk

Case Talk

Case Talk List


Case Talk Count : 2


Case Talk List

Status	Case Talk Name	Sharing Type	Date Shared	Expiration Date	
	test's Case_2020-09-20 10:50:46		2020-09-20	2020-12-19	
	test's Case_2020-09-20 10:49:27		2020-09-20	2020-12-19	



Go to Case Talk


Close


- Set share scope and access rights by clicking on the  icon.





Link Access
Who would you like to share the Case Talk link with?

 **Anyone** 
Anyone with a link can access

 **Share with a Medit Link Organization**
Any organization with a partnership or contact with My Organization

 **Designated person (patient consultation)**
The link will be shared via email. Please set a pincode to ensure the security.

Security Settings
Expiration date
 2020-12-19

 You can extend the period up to 90 days from today.

Modify

Link Access

- Anyone with a link
- Share with a Medit Link Organization
Any organization with a partnership or contract with My Organization
- Designated Person (patient consultation)
The link will be shared via email. Please set a pin code to ensure the security.

Security Settings

- Link expiration date
- Pincode settings

Other Settings

Choose if you'd like to

- Share the patient's name
- Allow comments


If you enable comments, you will be able to communicate through the shared link.

If you uncheck the option, the recipient can only view your comments but will not be able to add the ones of their own.

3.2.3 Web Viewer


Use a browser to view shared data in Case Talk.

It is possible to communicate with the users who have access to the link using annotations and comments.



► Annotation

You can add annotations to a specified location on the 3D model. The annotations you add can be viewed by anyone with the link.



- Click on “New Annotation” on an any specific spot.

Create Annotation

Please enter your annotation.


Please enter your message

Cancel

OK


Medit Link User Guide


- Type what you want to deliver, then click 'OK'. An annotation is generated.



- Enable the conversation through annotation.

Comments	Communicate with users who have access to the link.
User	Displays the currently connected users.
Settings	Change sharing settings.
Real-time screen sharing	Share your screen in real time.

- Click on the 'LIVE'  to share your Web Viewer screen with the connected users.




- Click **OK** to start sharing your screen.
- Other people who have accessed the link accept screen sharing.
- The screen of the user who requested screen sharing is synchronized with the screen of the person who received the request.

3.2.4 Model Tree

Provides the ability to group and display data according to attributes, view/hide data, adjust transparency, and so on.

Show/Hide

Click on the  icon to adjust the display status of your data.

Opacity adjustment

When you mouse over an element, a slider bar is displayed to adjust the transparency.

Use the slider bar to adjust the transparency of your data.








Extended menu

When you mouse over an element, the Extensions menu is displayed.



- Show This Only: Hides all other but the current one.
- Hide: Hides the current data.

▷ 3D Model Control

	Pan	Moves the model.
	Rotate	Rotates the model.
	Zoom in/out	Zooms in/out the model.
	Zoom Fit	Fits the model to the screen.
	Change Color Display Mode	Changes the color display mode for the model.

3.2.4.1 Recommended browsers

PC	Chrome, Firefox, Safari
Android	Chrome, Firefox
iOS	Safari, Chrome, Firefox

3.2.4.2 Minimum specifications for mobile devices

Device	Android Same or above the specs of Galaxy S7	iOS Same or above the specs of iPhone 7
OS	Android 7.1 Nougat	iOS 11.0
RAM	Above 3GB	Above 2GB
AP	Above Snapdragon 710	A10 Fusion


3.3 Partners

You can search for a partner to request a contract, and see which partners are currently contracted.

3.3.1 Partner Search


You can search for the Medit Link users based on their business name and/or address.

In the search results, you can request Partnership from an organization you'd like to work with.



3.3.2 Recommended Partners

Medit Link displays its users located close to you based on your address information.




3.3.3 My Partner List

My Partner List displays a list of the companies that you currently have a contract with.



On Partnership vs. Contract


- If you are in a partnership with a business, you can request a contract with the organization.
- Once you request a contract, you can proceed with payment for the placed orders.



3.3.4 Partner Invitation

You can proceed with a partnership request with a subscription request to a user who is not a member of Medit Link.

- Click '**Invite Partner**' to display the following screen:



Invite Partner

Invite your desired partner to join Medit Link. When the lab accepts your invitation, a partnership is automatically formed.

CancelSend Invitation

- Enter the email address of the vendor you want to request a partnership from and send the invitation mail to proceed with the invitation.
- If the invited company confirms the membership and proceeds with the membership, the partnership will be established, and the Medit Link service based on this will be available.

3.4 Membership

Medit Link allows to organize and share files via Cloud Storage.

From now on, the cloud storage capacity is provided depending on the type of your membership.

The available membership plans are

- Standard,
- Premium and
- Premium Plus.

Detailed information for each membership plan will be provided later, before the membership service is launched in the second half of 2021.

Before the launch of the service, all users are entitled to Premium Plus membership plan for free.

3.5 Redeem Code

Register the redeem code provided as a part of the promotion.

3.6 Payment


3.6.1 Payment Details (Clinics, Labs accounts)

You can check all payment history and transactions through the Payment Detail page.

Payment Detail						
Check the details of payment information here.						
Total Number of Transactions : 0 Case			2020-08-12 ~ 2020-09-10		All Partners	
Organization	Patient	Payment Method	Price	Payment Status	Transaction Date	Order ID
No results were found for your search.						

Medit Link User Guide

3.6.2 Payment account information (Lab account only)



This feature is only available to the administrator of Lab account.

- You must enter your payment account information so that your account can use the online payment system when ordering.




If you select a country that can support Stripe, you will be notified of the Stripe connection, and you can create a Stripe account by pressing the appropriate button.

3.7 My Information


3.7.1 My Account

You can view/modify the information you entered during the registration and check your membership.




3.7.2 My Organization (Admin account only)

View/modify information of your organization.



3.7.3 My Price List (Lab: Admin and Member accounts)

- Enter the products offered by your lab and the price information for each product.



Name	Type	Method	Material
Inlay Zirconia	Inlay	-	Zirconia
Inlay Gold	Inlay	-	Gold
Inlay Press Ceramic	Inlay	-	Press Ceramic
Inlay Resin	Inlay	-	Resin
Inlay Hybrid Ceramic	Inlay	-	Hybrid Ceramic

Medit Link **User Guide**



You can make a transaction based on the product list and price information if you have a contract relationship with the client.



Medit Link does not support price-based transactions for partnerships.

3.8 Help Centre

Once you press '**Help Center**' you get transferred to Medit's Help Centre which contains information about the company's products and their functions.

3.9 Announcements

Displays notifications for new products.

3.10 Notices and Updates

Displays notifications about using and updating the Medit Link service.


3.11 Download the App

- Download the installer for the latest version of the program.
 - You can install the program on another computer if you have it already installed after the purchase of the product.



Program updates are automatic when the Internet is connected.

4 Medit Link App Overview



1. Title Bar
2. Boxes
3. Detailed Functions
4. Case List

4.1 Title Bar

The title bar displays network status, cloud synchronization, membership grade, and notifications.

4.1.1 Network Status

Medit Link enables you to use all the features of the service seamlessly when you are connected to a network.

The Title Bar displays the network status for the smooth use of the service.

If the network connection is not good, the following restrictions occur:

- Some functions, such as ordering, creating and sharing Case Talk, cloud synchronization are restricted.
- The overall performance of the service can be reduced due to the impact of the network environment.

In this case, you can enter the Offline Mode to perform basic tasks, and then use the features that require internet once you are connected.

4.1.2 Cloud Synchronization

Shows the status of the original data that is not synchronized to the cloud and provides the ability to manage synchronization manually.

Scan Raw File Upload Status

Shows the uploading status of raw scan files.
The raw file contains all the information scanned by the scan program.
If you upload it to the server, you can resume the scan anytime and anywhere on a different PC.
You can upload the raw file automatically by turning on the 'Upload Raw Scan File Automatically' option in the settings.

Files Being Uploaded

Cancel Uploading

Case Name	Size
-----------	------

Local Disk Files

Upload All

Case Name	Size
-----------	------


Close

4.1.3 User Information

User information displays basic member information such as profile image, organization name, email address, and membership name. It also provides logout option and home page shortcuts.

My Information

1



2

Organization Name

Wonhoonchoi_Clinic_Dev

Email Address

wonhoon.choi@medit.com

Name

Wonhoonchoi_Clinic

3

Close

Log Out


4

Edit >

- 1 Displays the profile image.
- 2 Displays the organization name, email address, and membership name.
- 3 Logs out of the currently connected account.
- 4 Logs in with the currently accessed account information and access the homepage.
Check/modify/manage registered information


4.1.4 Notifications

Provides notifications related to contracts/orders and work processes.



- 1 Shows only the unread messages.
- 2 Quickly navigate to the specific case or webpage by selecting the corresponding notification.

4.2 Boxes Overview



Medit Link User Guide

- Boxes available in Lab and Clinic accounts.

		Case Box	Order Box	In Box	Work Box	Trash Box	App Box
Clinic Account	Admin	O	O	X	X	O	O
	Member	O	X	X	X	O (Restoring cases)	O
Lab Account	Admin	X	X	O	O	O	O
	Member	X	X	X	O	O (Restoring cases)	O

4.3 Settings

General Settings: Provides features for overall system setup.

The screenshot shows the 'General' settings page of the Medit Link application. On the left is a sidebar with 'Settings' and 'General' (selected). The main content area is divided into sections: 'Save Folder' with a text input 'D:/LinkSave/' and a folder icon; 'Language Settings' with a dropdown menu set to 'English' and an 'Update' button; 'Teeth Labeling' with radio buttons for 'FDI' (selected) and 'US'; 'Notification' with a checkbox 'Allow Notifications' (checked); and 'Uploading Raw Scan Files' with a checkbox 'Upload Automatically' (checked). At the bottom right of the main area is an 'Apply' button. At the bottom right of the entire window is a 'Close' button.

Save Folder

Specify a folder to store the results generated from Medit Link and connected programs.

Languages Settings

Select the program language.
Improve your current language and update in real time when a new language is added.


Teeth Labeling Notifications

Select a labeling system.
Set whether to receive notifications for program updates.

Raw Scan File Upload

Determines whether to automatically upload raw scan data to the cloud.
It helps to back up the original data and perform seamless scanning when using Medit Link with the same account on another PC.
The original data cloud backup only applies to data acquired by Medit Scan for Clinics.

Storage Management: Provides features for viewing and managing storage usage.



Available Space

Provides features for viewing storage usage.

Disk Cleanup

You can free up disk capacity by deleting data within a set period of time.


- Use the calendar to proceed with the search.
- Select the type of file you want to delete.
- Click [Delete Files](#) to delete the selected file.

Cases in Trash Box

Free the Trash Box to free up your disk.
Data deleted using this feature cannot be recovered.

Medit Link User Guide


Set Connected Program: Set/manage programs that work with Medit Link.



You can connect the following programs:

	Clinic Account	Lab Account
Scanning S/W	Medit Scan for Clinics	Medit Scan for Clinics Medit Scan for Labs
CAD	exocad	exocad Maestro 3D

Export Option: Set how to export files from Medit Link.



CAD Type

Select the CAD program you want to use for your design. The program will adjust the files based on this selection, allowing you to proceed with the design without any additional axis alignment.

Export Files To

Specify the path where you want to store the results.

Mesh Format

Specify the format of the final mesh file.

Image Format

Specifies the format of the exported image files.

Compress & Combine

Compress: Compresses all results into a single file.
Combine (mesh): Merges the results of both arches into one result file. Useful when importing several files into a program is not possible.

Naming Rule

Define the rules for the output file name:

- o Case creation date/scan date/order date/patient name/case name/order number.

Case Converting Tool: Used to move cases between different accounts.

Settings

General

Storage Management

Set Connected Program

Export Options

Case Converting Tool

Proxy Option

About

User Guide

Export

2020-09-13 ~ 2020-09-20

All

Search / Contains

Files to Export

Case Name	Patient Name	Date Modified
<input type="checkbox"/> test's Case	test	2020-09-20 오전 10:42
<input type="checkbox"/> test's Case - Clone - Clone - Clone	test	2020-09-18 오후 7:00
<input type="checkbox"/> test's Case - Clone	test	2020-09-18 오후 6:58
<input type="checkbox"/> test's Case - Clone - Clone	test	2020-09-18 오후 6:58
<input type="checkbox"/> test's Case - Clone	test	2020-09-18 오후 6:58
<input type="checkbox"/> aaaa1111	test	2020-09-18 오후 6:37
<input type="checkbox"/> 배송분리시스템2	test	2020-09-18 오후 6:14
<input type="checkbox"/> complete complete	S-9_Scan_CAMOrder_Test...	2020-09-18 오후 3:57
<input type="checkbox"/> complete complete - Clone	S-9_Scan_CAMOrder_Test...	2020-09-18 오전 10:19

File Name

CaseBox_2020-9-20_123057

Export Files To

C:/Users/Medit-Wonhoonchoi/Documents/

Export

Import

A new case is created for every imported case. It contains case data, patient information and scan file, but not the order information.

Import from

C:/Users/Medit-Wonhoonchoi/Documents/

Import

File List

File Name	Import Time
CaseBox_2020-9-15_104846.meditex	2020-09-16 오전 7:39
CaseBox_2020-9-15_104846.meditex	2020-09-16 오후 12:09
CaseBox_2020-9-14_143157.meditex	2020-09-18 오전 8:51
CaseBox_2020-9-14_143157.meditex	2020-09-18 오전 8:55

Close

Proxy Options: In some limited network environments, you can change network settings.

Settings

General

Storage Management

Set Connected Program

Export Options

Case Converting Tool

Proxy Option

About

User Guide

Proxy Server

Use Proxy Server ☒

Proxy Server IP (or DNS)

Proxy Server Port

0

Proxy Server Scheme

☐ HTTP
 ☒ HTTPS

Proxy Server Username (Option)


Proxy Server Password (Option)

Please apply the changes.

Apply

Close

About: Displays the details of Medit Link.



User Guide: Gets the user guide for Medit Link. Requires an Internet connection.

4.3.1 Patient Management

Provides the ability to add/edit/delete and merge the existing patients.

4.3.2 Notices and Updates

Displays program updates and critical alerts.

4.3.3 My Account


Transfers you to Medit Link homepage under the account that you are currently logged in with.

4.3.4 Help Center

Goes to the Help Center page. You can get various information about using Medit Link.

4.4 Detailed Functions of Case Box/Work Box

Case Box/Work Box allows to create cases and filter them using various options.





4.4.1 Creating cases

Create a case by pressing the 'New Case' button in the top right corner and entering information such as patient and case name.

4.4.2 Change viewing options (Clinic Account only)

Change how you view cases between gallery and list forms.



- 1 Case name
- 2 Patient name and birth date
- 3 Form information
- 4 Date Scanned
- 5 Case status

Filtering Options of Case List

Filtering using status

Quickly view/hide cases of the selected status.

Filtering by groups

View cases by groups.

Filtering using period settings


Displays cases within a specified time period.

Filtering by search

Use a variety of search criteria to quickly locate a case.

5. Medit Link App – Box Detail

5.1 Case Box / Work Box










- 1 Displays the case name.
- 2 Displays the key available features: scanning, CAD, ordering, Apps.
- 3 Tabs for each of the tasks.
- 4 Functions for managing cases and files (case deletion, sharing, file attachment, export).
- 5 Selects the scan target – fill arches, or specific teeth.
- 6 Allows to select type/style/material/shade information for form based on the CAD settings.
- 7 Displays the type/style/material/shade information applied to the selected teeth.

About Tabs

Scan →	CAD →	File Viewer
Enter information for basic scan tasks (form and scanning options). The information is passed to the connected scan program and used as the order information.	Set CAD options. The information you set up is passed to the CAD program.	View result files generated by scanning programs and CAD. Provides the ability to effectively display and manage 3D and 2D data, both acquired and attached.

About detailed functions of the case

Case Talk 	Creates a shared link for the case. Users accessing shared links can easily view 3D data using a browser such as Chrome, Firefox, and Safari without installing additional programs. You can set the sharing scope, access rights, and so on for the links.
Export 	Exports scan/CAD result files in the specified file format. Exported files can be used for future processes such as CAD, CAM, etc.
Attach 	Attaches files in various formats.
Clone 	Clones a case to the Case Box. Useful for replicating treatment information for the same patient, or for additional orders after replicating all information/files in the same case.
Delete 	Deletes the case (only for the cases created in Online Mode). You cannot delete cases that have been ordered.
Support Request 	Asks the Medit Technical Support Center for technical support.
Case History 	Provides detailed information about the case.

History

Date Modified

2020-09-21 오후 5:29

Date Scanned

Date Created

2020-09-21 오후 3:02

Case Size

0MB

Close

Memo



Shows the notes for the case.

The memos will be included during ordering and can be used to communicate more detailed information about the case.

Go to Detail



Provides a shortcut for the same case in another box (Case Box ↔ Order Box; In Box ↔ Work Box).

Status

Displays the status of the case.

5.2 Order Box / In Box

Medit Link 2.3.0 Case (Medit Link 2.3.0)

Order | File Viewer

Order ID: 5718378

Patient Name: Medit Link 2.3.0 | Partner Name: WonhoonChoi_Lab_Stage

Date Scanned: 2020-09-14 오전 7:43 | Date Ordered: 2020-09-20 오후 1:11

Requested Delivery Date: 2020-09-27 | Expected Delivery Date: 2020-09-27

Memo

No.	Type	Method	Material	Shade	Price
12	Inlay	-	Zirconia	A3	0.00
13	Inlay	-	Zirconia	A3	0.00
14	Inlay	-	Zirconia	A3	0.00
15	Inlay	-	Zirconia	A3	0.00
16	Inlay	-	Zirconia	A3	0.00
45	Onlay	-	Zirconia	A3	0.00
46	Onlay	-	Zirconia	A3	0.00
47	Onlay	-	Zirconia	A3	0.00
48	Onlay	-	Zirconia	A3	0.00
23	Implant Crown	Screw Type	Zirconia	A3	0.00
24	Implant Crown	Screw Type	Zirconia	A3	0.00
25	Implant Crown	Screw Type	Zirconia	A3	0.00
26	Implant Crown	Screw Type	Zirconia	A3	0.00
Total Price					USD 0.00

- 1 Displays the case name and patient name.
- 2 Provides functions for order management (accepting/cancelling; delivering/receiving).

- 3 Tabs for each of the tasks.
- 4 Shows case status and shortcuts.
 - Go to Detail: A shortcut to check the details of ordered cases. Click the button to go to the Case Box or Work Box case.
 - Status: Displays the job status of the case.
- 5 Displays order number, patient name, account name, scan date, order date, memo, etc.
- 6 Displays the product information and delivery information.



About Tabs

Order →


Check the order number, patient name, account name, scan date, order date, memo, and other detailed information related to the order.

File Viewer

View files generated by scan/CAD operations. It effectively displays and manages 3D models and 2D images, both acquired or attached.

5.3 App Box

In addition to the built-in programs, we offer a variety of additional ones. Download the programs of your liking and use them freely.



5.4 Trash Box

Trash Box is a space for keeping deleted cases before they are completely removed or recovered.



Data deleted from the Trash Box is deleted from the PC and server and cannot be recovered.

The screenshot displays the 'Trash Box' interface in the Medit Link application. The interface includes a sidebar with navigation icons for Dashboard, Case Box, Order Box, App Box, and Trash Box. The main area shows a table of deleted cases. At the top, there are buttons for 'Restore' and 'Delete Items', and a 'Selected Count 3' indicator. A search bar is also present. The table has the following columns: Status, Case Name, Patient Name, Lab Name, Deleted Date, Deleted by, and Size. Three items are selected for deletion, indicated by blue checkmarks in the Status column.

Status	Case Name	Patient Name	Lab Name	Deleted Date	Deleted by	Size
<input checked="" type="checkbox"/>	JHTEST's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	0
<input checked="" type="checkbox"/>	test's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	1417711732
<input checked="" type="checkbox"/>	test's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	0
<input type="checkbox"/>	tes's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	0
<input type="checkbox"/>	Test-Creation's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	16243482
<input type="checkbox"/>	test's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	0
<input type="checkbox"/>	Test-Creation's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	0
<input type="checkbox"/>	test's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	0
<input type="checkbox"/>	test's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	13728200
<input type="checkbox"/>	test's Case	test		2020-09-01 오전 8:59	Wonhoonchoi_Clinic	0


6. Detailed Explanation of the Process

6.1 Managing cases: creating, editing, cloning, deleting

6.1.1 Creating cases

▷ Creating cases in Case Box (Clinic)

- Click 'New Case' button in the upper right corner.



New Case Information


Patient Name*

Case Name*

Creating patient information involves consent from the patient.
Patient information can be shared with 3-parties on order.

[Terms & Conditions](#)

Cancel Register Register & Scan

- Enter basic information such as patient name.
 - Entering a name that is not in the patient list creates a new patient. You can enter a name in the patient list to select from the list.
- In case you're registering a case for an already existing patients, click the  button.
 - You can select a new patient after registering or editing information for an existing patient.

Patient List

Select a patient in the list to set the patient of the case. You can add or manage the patient information in Medit Link.

All Search / Contains + New Patient


	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Registered Date	Recent Visit Date	
<input type="radio"/>	test				0	2020-08-28	2020-08-28	
<input type="radio"/>	test				0	2020-08-28	2020-08-28	
<input type="radio"/>	test				13	2020-08-23	2020-08-28	
<input checked="" type="radio"/>	test				85	2020-08-23	2020-05-08	
<input type="radio"/>	Wonhoon Choi				1	2020-04-01	2020-04-01	
<input type="radio"/>	ttttt				0	2020-08-18		
<input type="radio"/>	13				0	2020-08-20		
<input type="radio"/>	14				0	2020-08-20		
<input type="radio"/>	15				0	2020-08-20		

Close Select

- Case name will be filled in automatically, however you can still change it.
- Press 'Register & Scan' to complete the case registration and proceed to scanning immediately.
 - Otherwise, press 'Register' which is useful when registering multiple cases at the same time.
 - You will be redirected to the main page of the Case Box.

▷ **Creating cases in Work Box (Lab)**

- Click '**New Case**' button in the upper right corner.



The screenshot shows a form titled "New Case Information". It contains four input fields: "Clinic Name*" (a dropdown menu), "Patient Name*" (a text box), "Case Name*" (a text box), and "Desired Delivery Date" (a date picker showing "2020-09-08"). At the bottom right of the form are two buttons: "Cancel" and "Register".

- Enter or select the clinic name.
- Enter the patient name and case name.
 - Case name will be filled in automatically, however you can still change it.
- Enter the expected delivery date considering the completion date of the work.
- After registering the case, go to the job information entry screen.


6.1.2 Editing cases

▷ Editing cases in Case Box (Clinic)

- Click on the  icon next to the case name.


Edit Case Information

Patient Name*




✓ 'test' is selected.

Case Name*

 Creating patient information involves consent from the patient.
 Patient information can be shared with 3-parties while placing an order.

[Terms & Conditions](#)

Cancel
Apply










- Click on the  icon to change the patient name. Note that it cannot be edited; instead, you can find a patient name you'd like to replace the current one with.

Patient List


Select a patient in the list to set the patient of the case. You can add or manage the patient information in Medit Link.

All

+ New Patient


	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Registered Date	Recent Visit Date	
<input type="radio"/>	test				0	2020-08-28	2020-08-28	
<input type="radio"/>	test				0	2020-08-28	2020-08-28	
<input checked="" type="radio"/>	test				13	2020-08-23	2020-08-28	
<input type="radio"/>	test				85	2020-08-23	2020-05-08	
<input type="radio"/>	Wonhoon Choi				1	2020-04-01	2020-04-01	
<input type="radio"/>	ttttt				0	2020-08-18		
<input checked="" type="radio"/>	13				0	2020-08-20		
<input type="radio"/>	14				0	2020-08-20		
<input type="radio"/>	15				0	2020-08-20		

Close
Select

- Check the patient name you'd like to replace the current one with.
- To create a new patient, press the 'New Patient' button.
- Press  to edit the registered information about the patient.
- Press 'Apply' to finish the process.

▷ Editing cases in Work Box (Lab)

- Click on the  icon next to the case name.



The form is titled "Edit Case Information". It contains four input fields: "Clinic Name*" with a dropdown menu showing "test", "Patient Name*" with a text input showing "test", "Case Name*" with a text input showing "test's Case", and "Requested Delivery Date" with a date picker showing "2020-09-22". At the bottom right, there are two buttons: "Cancel" and "Apply".


- Press 'Apply' to after the editing is done.


 You cannot edit the information of the cases received through the order in In Box.

6.1.3 Cloning cases

You can use the Clone feature to clone a case to rescan with the same treatment information. You can duplicate the entire data to order a case that has already been ordered, or chose other cloning options.

▷ Cloning cases in Case Box (Clinic)

- Click on the  icon.
- You will be asked what kind of information you'd like to clone to the new case.




The dialog box is titled "Information". It contains an information icon (i) and the text "Select a clone type.". At the bottom, there are three buttons: "Cancel", "Form Data Only", and "All Data".

- Choose whether to replicate only treatment information or data all together.


- When cloning is complete, a case with the name "[Case name] – Clone" will be created.


▷ Cloning cases in Work Box (Lab)

- Click on the  icon.
- Duplicate all data to create a new case.
- When cloning is complete, a case with the name "[Case name] – Clone" will be created.

6.1.4 Deleting cases (Lab and Clinic)

You can delete cases.

- Click on the  icon.
- Click 'OK' to delete the case.



- Deleted cases are moved to the Trash Box and can be recovered if necessary.

6.2 Inputting form information

You can enter treatment information based on the connected program of your selection.
The information you enter will be used as basic information for scanning, CAD and ordering.

6.2.1 Inputting form information in Case Box

▷ Entering dental information in the Scan tab

- Select whether you are going to input treatment information by teeth or by arch.
- Select Type – Method – Material – Shade.
 - If you entered information incorrectly, you go back to select it again.
 - The selected information is automatically registered to the order list.
- Click '**Save**' to save your registration.

▷ Entering dental information in the CAD tab

- Select one of the following options that will define your scanning strategy: Teeth/Flexible Multi-die/Multi-die.
- Select the tooth number/arch/die number.
- Define the settings for detailed options and CAD operations.
- Click '**OK**' to save the settings.
 - The selected information is automatically registered to the order list.
- Click '**Save**' to save your registration.

6.2.2 Inputting form information in Work Box

▷ Entering dental information on the Scan tab

- Select the scan program you want to use.
- Select one of the following options that will define your scanning strategy: Teeth/Flexible Multi-die/Multi-die.
- Methods may differ depending on the selected scan program.

	Medit Scan for Clinics	Medit Scan for Labs
Teeth	<input type="radio"/>	<input type="radio"/>
Arch	<input type="radio"/>	<input type="radio"/>
Flexible Multi-die	<input checked="" type="radio"/>	<input type="radio"/>
Multi-die	<input checked="" type="radio"/>	<input type="radio"/>

- Select teeth, arch or die numbers for the scanning process.
- Select Type – Method – Material – Shade.
 - If you entered information incorrectly, you go back to select it again.
 - The selected information is automatically registered to the list.
- Click '**Save**' to save your registration.

▷ Entering dental information in the CAD tab


- Select one of the following options that will define your scanning strategy: Teeth/Flexible Multi-die/Multi-die.
- Select the tooth number/arch/die number.
- Define the settings for detailed options and CAD operations.
- Click '**OK**' to save the settings.
 - The selected information is automatically registered to the list.
- Click '**Save**'.

6.2.3 Inputting bridge information


You can enter information for a bridge in case same information is applied to all selected teeth.

- When selecting several teeth of the same type a bridge icon will appear by their side. You can turn the bridge on/off by clicking the icon.

Before setting a bridge



After setting a bridge



6.2.4 Setting scanning options

Depending on the program you are working with, various scan options are provided:
The selected scan options are passed to each scan program.

▷ Medit Scan for Clinics



Pre-op for Maxilla

When checking this option, you will be able to scan maxillary diagnostic (pre-operation) model.



Pre-op for Mandible

When checking this option, you will be able to scan mandibular diagnostic (pre-operation) model.



Face Scan

When checking this option, you will be able to acquire face scan data such as teeth, lips, and nose with the intraoral scanner.

▷ Medit Scan for Labs

➤ General options



Adjacent teeth

Mark the prepared and adjacent teeth and pass this information to the scanning program.



Antagonist

Automatically marks the antagonist and passes this information to the scanning program.



Articulator

Displays options for selecting the articulator in the scanning program.



Movement
Market

Allows to use the Jaw Motion scan strategy in the scanning program.

➤ Options by teeth



Scanbody

Allows to set the strategy for scanning scanbody in the scanning program.



Pre-op Scan

Allows to set the strategy for scanning pre-op model in the scanning program.



Wax-up
Scan

Allows to set the strategy for scanning wax-ups in the scanning program.


6.3 Copying form information

You can copy the treatment information.

▷ Copy dental information using context menu

- Select the tooth for which the tooth information is set.
- Right-click and select '**Copy**' from the context menu.


Medit Link User Guide




- Select the target tooth.
- Click  to complete the copy.

▷ Copy dental information using shortcut keys


- Select the tooth for which the tooth information is set.
- Press and hold CTRL and select the target tooth.

6.4 Deleting treatment information

- You can delete the dental information by clicking  in the order list.

Order List					
No.	Type	Method	Material	Shade	
16	Inlay	-	Zirconia	A4	
25	Inlay	-	Gold	-	
26	Onlay	-	Gold	-	

▷ Delete dental information using context menu

- Select the tooth for which the tooth information is set.
- Right-click to select "Delete" from the context menu.
- Select the target tooth.
- Click  to complete the copy.

▷ Delete dental information using shortcut keys


- Select the tooth for which the tooth information is set.
- You can delete the selected information by pressing the keyboard "Del" button.

6.5 Adding product information to 'Favorite'

You can quickly select a product by registering it as a favorite.

The list of registered favorites is managed as a template and can be shared through the export and import features.

6.5.1 Adding to favorites

- Select the tooth number.
- Select Type – Method – Material – Shade.
 - Type-style-Material-Shades can be registered as favorite even if not all of them are selected.
- Click  .
- A combination of the selected parameters, excluding tooth number, is listed on your favorite search list.

6.5.2 Managing favorite templates

You can register, manage, and share different types of templates.


Add	Adds a new template.
Rename	Renames the current template.
Delete	Deletes the current template.
Export	Exports the current template to a file.
Import	Imports a template file.

6.6 Copying form information between Scan and CAD tabs

You can copy the treatment information entered between Scan tab and CAD tab.

6.6.1 From Scan tab to CAD tab

- If there is no information input in the CAD Tab, it automatically copies the information when you go from the Scan Tab to the CAD Tab.
- If there is information in the CAD tab, you will have to copy it manually.
 - Right-click on the CAD tab.
 - Select **"Convert form Information from Scan"** from the context menu.




- The information in the Scan tab is copied to the CAD tab.

6.6.2 From CAD tab to Scan tab

To copy information, use the following methods:

- Right-click on the Scan tab.

- Select "**Convert form Information from CAD**" from the context menu.



 Note that only Type-Method information will be copied.

6.7 Scanning process

Based on the entered information, you proceed with scanning.

- Enter/save treatment information and select the necessary scan options.
- Click '**Scan**' button in the upper-right corner of the window.

6.8 CAD process

Medit Link is integrated with exocad and Maestro 3D CAD programs.


CAD	Registering treatment information	Confirmation of treatment information	Direct CAD Execution	CAD Results import back to Medit Link
exocad	O	O	O	O
Maestro 3D	X	X	O	X

- In Settings → Set Connected Program, you can link the CAD program of your choice from the list.
- Enter/save task information and design options in the CAD Tab.
- Click '**CAD**' in the upper-right corner of the window.
- The linked CAD program is run and the data is sent automatically.
- Save the results after performing CAD.

6.9 Exporting result files / Case Talk

6.9.1 Export

In addition to sharing results with partners through orders, Medit Link allows you to store your results in the desired file format using the export function.



Export Options

- Set Export Folder Name**
2020-09-18-ttt
- CAD Type**
exocad
- Export Files To**
C:/Users/Medit-Wonhoonchoi/Desktop/
- Mesh Format**
☐ Keep Original Mesh Format
☒ Mesh Format
☒ OBJ ☐ PLY ☐ STL
- Image Format**
☐ Keep Original Image Format
☒ Image Format
☒ JPG ☐ JPEG ☐ PNG ☐ BMP
- Combine**
☐ Combine Individual Mesh
- Files to Export**

Name	Type	Size
<input checked="" type="checkbox"/> 5657377.colLabProject	colLabProj...	3 KB
<input checked="" type="checkbox"/> 5657449.dentalProject	dentalProj...	33 KB
<input checked="" type="checkbox"/> 5657449-13	meditMesh	258 KB
<input checked="" type="checkbox"/> 5657449-14	meditMesh	279 KB
<input checked="" type="checkbox"/> 5657449-15	meditMesh	277 KB
<input checked="" type="checkbox"/> 5657449-16	meditMesh	271 KB
- Compress Files**
☐ Compress Files

Close Export

Set Export Folder Name

Specify the folder name. The exported file will be named according to the naming rule used by the selected CAD program.

CAD Type

Select the CAD program to be used for your design.

Medit Link will adjust the files based on this selection, allowing to proceed with the design without any additional axis alignment.

Export Files To

Specifies the path where you want to store the results.

Mesh Format

Specifies the format of the mesh file.

Image Format

Specifies the format of the image file.

Combine

Merges the results of both arches into one result file.

	Useful when importing several files into a program is not possible.
Files To Export	Displays and allows to set the list of files to export.
Compress Files	Check this option to compress the results to a single file.

6.9.2 Case Talk

You can share your work with multiple people using Case Talk and view shared files through Internet browsers such as Chrome, Firefox, and Safari without installing a separate program.

Share to Case Talk

[New Case Talk](#) [Case Talk List](#) Case Talk Count : 0

File List

	Name	Type	Size
<input checked="" type="checkbox"/>	Maxilla	meditMesh	763 KB
<input checked="" type="checkbox"/>	Default_1.stl	stl	2,694 KB
<input checked="" type="checkbox"/>	Default_2.stl	stl	2,596 KB
<input type="checkbox"/>	Default_3.stl	stl	2,593 KB

Get Link

3 items selected

5.91MB / 200MB

If you share files over 200MB, Web Viewer may not work on some mobile devices.
Please check the recommend requirement

[Check requirement](#)

You can create Case Talk item and share the link with anyone, specific people, organization in Medit Link.
You can set the password through website to limit the access to shared link.

[Go to Case Talk](#)

Please note that Medit Link is not responsible for any problems caused by sharing the public link you created.

[Close](#)

► Case Talk

- Select the files to share.
- Click '**Get Link**' to create a shared link.
- The created shared link can be found in the Case Talk list.

Case Talk Toolbox



Opens the Case Talk directly in the Web Viewer.



Deletes shared items.



Change sharing settings.

You can set the access rights, and so on.

6.10 File Viewer

File Viewer groups and shows various data, including the data acquired from programs linked to MEDIT Link and data attached by users.

It provides basic management functions such as viewing/hiding data, adjusting transparency, and deleting data.

It supports previews of files in the following formats:

- 3D data: STL, OBJ, PLY
- 2D images: PNG, JPG, JPEG, BMP, GIF, TIF, TIFF
- Video files: MKV, MP4, AVI, WMV

6.11 Model Tree

It provides the ability to group and display data and perform actions such as viewing/hiding, transparency adjustment, and deletion through the model tree.

Show/Hide

You can adjust the view/hide status of your data.



Transparency
adjustment


You can adjust transparency for elements in the "View" state by moving the slider bar.



Extension Menu

Displays additional features available.





- **Show this only:** Hides all other data and display only current data.
- **Hide:** Changes the current data to the Hide state.
- **Zoom fit:** Zooms the screen around the current data.
- **Delete:** Deletes the data.

▷ 3D Model Control


	Pan	Moves the model.
	Rotate	Rotates the model.
	Zoom in/out	Zooms in/out the model.
	Zoom Fit	Fits the model to the screen.
	Change Color Display Mode	Changes the color display mode for the model.

▷ 2D Images Control

	Next/previous image	Moves among images.
	Edit	Edits the image.
	Pinch/Paste Window	You can detach the image window and adjust it to the desired position. You can resize the window. When the window is restored to its original location, it returns to its default size.

6.11.1 Covering empty spaces


Creates a model by covering the empty spaces.



6.11.2 Capture and Memo


Capture the current screen as a 2D image, draw on it, or take notes.

You can save it as a file together with your notes and save them as a reference.



6.11.3 Orthodontic Capture

Captures the images in current viewing direction and top/bottom/left/right from the current viewing direction.



6.12 Order process and order management

You can easily put both acquired and attached data into an order to your partner.

6.12.1 Ordering

You can order cases to your partners from Case Box.

- Click '**Order**' button.
- Check the information required for your order.

Order Case


Partner*

WonhoonChoi_Lab_Stage

Go To Partner

Payment Details

File For Order



No Price Information

Case Name*

test's Case

Patient Name*

test

☒ Share Patient Name

Delivery Desired Day

2020-09-09

오전 12:00

Memo

Cancel

OK


- Select the partner's account.
 - Press '**Go to Partner**' to view information about selected partners.
- Check the product information and files.
 - You can decide whether to share patient name.
- Enter the requested delivery date.
 - The default value is set to 7 days after the order date.
- You can memos to your partner to outline some specific details.
- Proceed with your order.
 - **Skip payment**: Proceed with the order without using the payment system provided by Medit Link. Available when ordering from a partner with the 'Contracted' status.
 - **Pay**: Utilize the payment system provided by Medit Link to proceed with the order. Available when ordering from a contracted partner.
 - **OK**: We proceed with the order without any conditions regarding the price of the product. Available when ordering from a partner with the 'Partnership' status.

Medit Link **User Guide**

▷ Confirm Order – In Box

➤ Confirm Order


- Accept the order upon receiving it.
- Accepted orders will be forwarded to the Work Box.



▷ Cancel/Reject Order – Order Box/In Box


➤ Cancelling Order – Order Box

- You can cancel an order before your partner accepts it.
- Canceled cases can be re-ordered.




➤ Reject order – In Box

- You can reject a received order.
- Rejected cases can be re-ordered.



6.12.2 Job Completed – Work Box

- Once all tasks are done, press "Done" to complete the case work.




➤ Enter shipping information – In Box

- Enter the shipping information for the completed task and begin shipping.

➤ Enter shipping information

- Click on the "Shipment" tab.



test's Case / Wonhoonchoi-Clinic

Order ID: 7991903

Patient Name: test Partner Name: Wonhoonchoi-Clinic

Scanned Date: 2020-09-07 오후 2:08 Ordered Date: 2020-09-09

Delivery Desired Day: 2020-09-14 12:00 오전 Delivery Scheduled Day: 2020-09-09

Memo: <Scanbody Library>
 2 : DAS_PREMIUM_L12_0003 - TI-BASE NR (GH=0.3) (CH=5) - 0°
 3 : DAS_PREMIUM_L12_0003 - TI-BASE NR (GH=0.3) (CH=5) - 0°
 4 : DAS_PREMIUM_L12_0003 - TI-BASE NR (GH=0.3) (CH=5) - 0°
 5 : DAS_PREMIUM_L12_0003 - TI-BASE NR (GH=0.3) (CH=5) - 0°
 6 : DAS_PREMIUM_L12_0003 - TI-BASE NR (GH=0.3) (CH=5) - 0°

No.	Type	Method	Material	Shade	Price
-1	Type	-	-	-	■

Total Price: KRW 500

Medit Link User Guide

- Enter shipping information:

Shipment Registration

Select a delivery method and enter the required information to register.
The information you enter is saved in the Shipment tab and can be checked at any time.

☒ Shipment ☐ Self-Delivery

Company

Tracking Number

Memo

Cancel

Register

- Click the "**Register**" button to complete the shipment information entry.
 - The shipping information you entered will be forwarded to your partner.
- Confirmation and Completion of Shipping – Order Box
 - After receiving the product proceed with the delivery confirmation.
 - The status of the case changes to Completed and the entire order process is completed.

<

test's Case - clone (test)

Confirm Delivery


Order

File Viewer

Shipped

6.12.3 Checking order information

You can check the details of the case you ordered from Order Box in Clinic Account and in the Order Tab in your Lab Account.



The screenshot shows the Medit Link interface with a sidebar on the left containing icons for Dashboard, Case Box, Order Box, App Box, and Trash Box. The main area is titled 'test's Case (test)' and has two tabs: 'Order' and 'Result'. The 'Order' tab is active and contains a form with the following fields:

- Order ID: 3705170 (marked with a blue circle 1)
- Patient Name: test
- Partner Name: Wonhoonchoi_Lab_Dev
- Scanned Date: (empty)
- Ordered Date: 2020-08-11 오전 8:39
- Delivery Desired Day: 2020-08-18, 12:00오전
- Delivery Scheduled Day: 2020-08-17
- Memo: (empty text area)

On the right side of the 'Order' tab, there is a 'Shipment' section (marked with a blue circle 2) containing a table with the following data:

No.	Type	Method	Material	Shade	Price
23	Inlay	-	Zirconia	A3	0.00
24	Inlay	-	Zirconia	A3	0.00
25	Inlay	-	Zirconia	A3	0.00
26	Inlay	-	Zirconia	A3	0.00
27	Inlay	-	Zirconia	A3	0.00
28	Inlay	-	Zirconia	A3	0.00

At the bottom of the 'Shipment' section, the 'Total Price' is listed as USD 0.00.

- Order information
 - The information that can be included in the order is as following: Order ID/patient name/account name/scan date/order date/delivery date/delivery due date/memo
 - You can check the details of the ordered product and verifying the delivery information.
- Product list displays the details of the ordered product.
 - Shipping: Enter/check the shipping information.

6.12.4 Managing expected/requested delivery dates

You can enter the expected delivery date when placing an order, and this information will be forwarded to your partner.

- Enter requested delivery date – Clinic Account
 - When ordering from Case Box, you can enter the requested delivery date.

Medit Link User Guide

Order Case

Partner*
Wonhoonchoi_Lab_Dev [Go To Partner](#)

Payment Details **File For Order**

No Price Information

Case Name*
test's Case

Patient Name*
test

☒ Share Patient Name

Delivery Desired Day
2020-09-08 오후 12:00

Memo

[OK](#) [Close](#)

- The default value is 7 days after the date of the order.
 - Enter the date/time you wish to receive the products and proceed with the order.
- Enter expected delivery date – Lab Account
- You can enter the expected delivery date for the case ordered from In Box.

Medit Link Wonhoonchoi_Lab_Stage

Medit Link2.3.0 Release / wonhoonchoi Beta [Reject](#) [Confirm Order](#) Pending

Order **File Viewer**

Order ID
5718378

Patient Name
Medit Link 2.3.0

Partner Name
wonhoonchoi Beta

Date Scanned
2020-09-20 오후 1:11

Date Ordered
2020-09-20 오후 1:11


Requested Delivery Date
2020-09-27 오후 01:10

Expected Delivery Date
2020-09-27 오후 01:10

Memo

Order List **Shipment**

No.	Type	Method	Material	Shade	Price
12	Inlay	-	Zirconia	A3	0.00
13	Inlay	-	Zirconia	A3	0.00
14	Inlay	-	Zirconia	A3	0.00
15	Inlay	-	Zirconia	A3	0.00
16	Inlay	-	Zirconia	A3	0.00
45	Onlay	-	Zirconia	A3	0.00
46	Onlay	-	Zirconia	A3	0.00
47	Onlay	-	Zirconia	A3	0.00
48	Onlay	-	Zirconia	A3	0.00
23	Implant Crown	Screw Type	Zirconia	A3	0.00
24	Implant Crown	Screw Type	Zirconia	A3	0.00
25	Implant Crown	Screw Type	Zirconia	A3	0.00
26	Implant Crown	Screw Type	Zirconia	A3	0.00
Total Price					USD 0.00

- Check the requested delivery date with the partner and make adjustments using .
- Enter the expected delivery date considering the actual working day, then click 'Apply'.

< **September** **2020** >

S	M	T	W	T	F	S
30	31	1	2	3	4	5
6	7	8	9	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30	1	2	3
4	5	6	7	8	9	10

Delivery Time

오전 12:00
^
v


Comment

Test - Desired Delivery Day |

Cancel
Apply

The expected delivery date and memo will be forwarded together with the order.

- You can check the date and memos in the Order Box.
- You can change the expected delivery date, and the changes will be forwarded to your partner.

- You can check the case history using .

6.13 Patient Management

Patient management tool provides functions for adding/editing/merging patients.

Patient Management

Allows to select patients to delete or merge their information from different cases.
 Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
 Press 'Delete' to remove patients from the system.

All [+ New Patient](#)

Patients with No Cases ☐

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date	
<input checked="" type="checkbox"/>	5-9_JScan_CAMOrder_Teeth2				3	2020-09-18	2020-09-18	
<input checked="" type="checkbox"/>	5-9_JScan_CAMOrder_Teeth2				1	2020-09-18	2020-09-18	
<input checked="" type="checkbox"/>	test				2	2020-09-17	2020-09-17	
<input checked="" type="checkbox"/>	test				22	2019-07-26	2020-09-17	
<input checked="" type="checkbox"/>	test				1	2020-09-17	2020-09-17	
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16	
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16	
<input type="checkbox"/>	gltt				1	2020-09-15	2020-09-15	
<input type="checkbox"/>	Medit Link 2.3.0			2020-9-24	2	2020-09-14	2020-09-14	
<input type="checkbox"/>	test				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	테스트				2	2020-09-13	2020-09-13	
<input type="checkbox"/>	test3				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	test1				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	testq1231				26	2019-07-25	2020-09-11	

[Close](#)

6.13.1 Adding patients

- Click 'New Patient'.

New Patient

Patient Name*

General

☐ Female ☐ Male

Birth Date

Patient ID

Phone Number

Email Address

Memo

Creating patient information involves consent from the patient.
 Patient information can be shared with 3-parties on order.

[Terms & Conditions](#)

- Enter the required information in the pop-up window, then click 'Save' to save the information.

6.13.2 Deleting patients

- Select the patient you want to delete from the list.
 - Patient deletion is possible if there no cases to the patient.

Patient Management

Allows to select patients to delete or merge their information from different cases.
 Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
 Press 'Delete' to remove patients from the system.

All Search / Contains + New Patient

Delete Merge Patients Patients with No Cases

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date	
<input type="checkbox"/>	test				0	2019-08-28	2019-08-28	
<input type="checkbox"/>	test				0	2019-08-28	2019-08-28	
<input type="checkbox"/>	test				0	2019-08-28	2019-08-28	
<input type="checkbox"/>	A-I Abument2				0	2019-08-28	2019-08-28	
<input type="checkbox"/>	Scanbody				0	2019-08-28	2019-08-28	
<input type="checkbox"/>	Margin Line Creation				0	2019-08-28	2019-08-22	
<input type="checkbox"/>	Abument Library				0	2019-08-28	2019-08-22	
<input type="checkbox"/>	GuideMessage				0	2019-08-28	2019-08-16	
<input checked="" type="checkbox"/>	margin				0	2019-08-28	2019-08-07	
<input checked="" type="checkbox"/>	test2				0	2019-08-28	2019-08-07	
<input checked="" type="checkbox"/>	suljaskidf				0	2019-08-28	2019-08-07	
<input checked="" type="checkbox"/>	Test1				0	2019-08-28	2019-08-07	
<input checked="" type="checkbox"/>	test				0	2019-08-28	2019-07-28	
<input checked="" type="checkbox"/>	test				0	2019-08-28	2019-07-25	

Close

- Click '**Delete**'.

6.13.3 Merging patients

Manage all cases belonging to one patient all together by merging patient information.

- Select the patient you want to merge from the list.

Patient Management

Allows to select patients to delete or merge their information from different cases.
 Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
 Press 'Delete' to remove patients from the system.

All Search / Contains + New Patient

Delete Merge Patients Patients with No Cases

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date	
<input type="checkbox"/>	5-9 Jicam_CAMOrder_test2				3	2020-09-18	2020-09-18	
<input type="checkbox"/>	5-9 Jicam_CAMOrder_test2				1	2020-09-18	2020-09-18	
<input checked="" type="checkbox"/>	test				2	2020-09-17	2020-09-17	
<input checked="" type="checkbox"/>	test				22	2019-07-28	2020-09-17	
<input checked="" type="checkbox"/>	test				1	2020-09-17	2020-09-17	
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16	
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16	
<input checked="" type="checkbox"/>	gmt				1	2020-09-15	2020-09-15	
<input type="checkbox"/>	Medit Link 7.3.0			2020-9-24	2	2020-09-14	2020-09-14	
<input type="checkbox"/>	test				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	测试				2	2020-09-13	2020-09-13	
<input type="checkbox"/>	test3				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	test1				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	test1221				26	2019-07-25	2020-09-11	

Close

- Click '**Merge Patients**'.

Medit Link User Guide

- Select the patient you'd like to merge all cases to.
- All cases are merged based on the selected information.

Merge all patient information to

Please select the patient you'd like to merge all data to.
After merging all data together, all information in other patient files will be deleted.

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date
<input checked="" type="radio"/>	5-9_iScan_CAMOrd...				1	2020-09-18	2020-09-18
<input type="radio"/>	test				2	2020-09-17	2020-09-17
<input type="radio"/>	test				22	2019-07-26	2020-09-17

[Cancel](#) [Apply](#)

- Click '**Apply**' to proceed with the merge.

6.14 Case filtering

Use various filtering options to quickly identify the desired patient and case.

6.14.1 Filtering by status

Provides the ability to filter cases according to their status.

☒ All | ☐ Form | ☐ Scanned | ☐ CAD | ☐ CAM | ☐ Mill | ☐ Done

- You can click the status you want to view to display only cases in that state.
- You can select multiple states at a time.


Different filtering options for each box:


Case Box	Form/Scanned/CAD/Ordered/Completed
Order Box	Pending/Accepted/Shipped/Completed
In Box	Pending/Accepted/Ready/Shipped/Completed
Work Box	From/Scanned/CAD/Done
Trash Box	Not supported

6.14.2 Filtering using period settings

Filters and displays cases within a set period of time.

► To set the duration:

- Click  2020-08-10 ~ 2020-09-09.
- Select the date range you want to view.




- Custom Range: Specify the start and end dates to display cases within this period.
- Today: Displays the case for today's date.
- Yesterday: Displays the case of yesterday's date.
- Past 7 days: Shows cases before 7 days except today.
- Past 30 days: Shows cases before 30 days except today.
- This Month: Show this month's case.

The case corresponding to the date you changed is updated in real time.

6.14.3 Grouping cases

Use various grouping and sorting methods for your cases.



- Grouping Method: Displays the entire case grouped by the selected method.
- Sort by: Displays cases within a group sorted by the selected method.

For each box, the following grouping and sorting methods can be applied:

Box	Grouping	Sorting
Case Box	None	Patient Name / Date Modified / Date Scanned / Date Ordered
	Patient Name	Date Modified / Date Scanned / Date Ordered
	Date Scanned	Patient Name / Date Modified / Date Ordered
	Date Ordered	Patient Name / Date Modified / Date Scanned
Order Box	None	Patient Name / Date Scanned / Date Ordered/ Requested Delivery Date
	Date Scanned	Patient Name / Date Ordered / Requested Delivery Date
	Date Ordered	Patient Name / Date Scanned / Requested Delivery Date
	Requested Delivery Date	Patient Name / Date Scanned / Date Ordered
In Box	None	Patient Name / Date Ordered / Requested Delivery Date
	Date Ordered	Patient Name / Date Scanned / Requested Delivery Date

	Requested Delivery Date	Patient Name/ Date Ordered
Work Box	None	Patient Name / Date Modified / Requested Delivery Date
	Requested Delivery Date	Patient Name / Date Modified

 The Trash Box does not provide a collective view and sorting function.

6.14.4 Filtering with search capabilities

Leverage search capabilities to quickly display the information you want.

All

▼

Search / Contains

🔍

Each box provides the following search options:

Case Box	All, Case Name, Patient Name, Patient ID, Form Information, Lab Name / Contains, Equal
Order Box	All, Case Name, Patient Name, Patient ID, Lab Name, Order ID / Contains, Equal
In Box	All, Case Name, Patient Name, Clinic Name, Order ID / Contains, Equal
Work Box	All, Case Name, Patient Name, Form Information, Clinic Name, Order ID / Contains, Equal
Trash Box – Clinic	All, Case Name, Patient Name, Lab Name, Deleted By / Contains, Equal
Trash Box – Lab	All, Case Name, Patient Name, Clinic Name, Deleted By / Contains, Equal