



jota product line SINTERED DIAMOND ASSORTMENT

Improved Quality and Performance

- › Outstanding grinding performance and durability
- › Modified shapes and dimensions for precise results
- › No discoloration on ceramics

**valid from
13.10.2014**



M801			spherical (round)	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M801	HP	807 104 001 524 -		033
			L mm	3,3
				2

M807			inverted conical	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M807	HP	807 104 030 524 - 		025
			L mm	3,7
				2


M856			bullet, conical pointed, slender	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M856	HP	807 104 274 524 -		060
M856G	HP	807 104 274 534 -		060
			L mm	14,0
				2


M862			flame
			
Fig	Shank	ISO	Ø
Handstück Straight Handpiece			
M862	HP	807 104 250 524 -	037
		L mm	11,0
			2

M813			inverted conical	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M813	HP	807 104 014 524 -		025
			L mm	1,0
				2

M837			side and end cutting	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M837G	HP	807 104 113 534 -		050
			L mm	12,0
				2

M863			cylindrical, long ogival end	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M863	HP	807 104 250 524 -		016
			L mm	8,0
				2

M941			disc thin	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M941	HP	807 104 345 524 -		220
			L mm	0,3
				1

M852		conical domed	
			
Fig	Shank	ISO	Ø
Handstück Straight Handpiece			
M852	HP	807 104 167 524 -	035 050
		L mm	12,0 12,0
			2 2

M854			conical	
				
Fig	Shank	ISO		Ø
Handstück Straight Handpiece				
M854	HP	807 104 173 524 -		031
			L mm	10,0
				2

Read carefully before use:

- Clean before first use on ceramics with JOTA cleaning stone 529D
- Clean after use on metal alloys with JOTA cleaning stone 529D
- Periodical cleaning with JOTA cleaning stone 529D grants best cutting performance
- Do not exceed maximum speed of 25'000 rpm
- Do not apply excessive pressure during preparations to avoid breakage risk

529D			
			
Fig	Shank	ISO	
529D		635 000 600 524 750	75 x 25 x 8 mm
			
			1